

ISSUE 5 August 2010

www.afootprint.eu

Personalised biomechanical models for intelligent orthotic design and optimisation

The A-FOOTPRINT project team provides some early insights on the development of personalised biomechanical models of the lower limb and foot

Two work packages in the A-FOOTPRINT project are currently developing patientspecific biomechanical models of lower limb and foot. These models will be used in the optimisation design and testina of personalised foot and ankle-foot orthoses. Researchers from Glasgow Caledonian University (GCU) and Maastricht University Medical Center (MAS) are recruiting normal healthy adult subjects and patients with a range of foot and ankle pathologies to undergo a series of foot examinations, anthropometric measurements, computed





Fig 1. 3D Surface scan of the foot and lower leg




AnyBody Technology (ABT) and **MAS** are currently undertaking the musculoskeletal and orthotic modeling using the information described above. The models start as generic models,




A-FOOTPRINT is a thirteen partner Consortium comprising six SMEs, five HEIs, one large-industry, and one research organisation partners from seven EU member states. The project is funded under **FP7 Cooperation Work Programme: Theme 4 - Nanosciences**, **Nanotechnologies, Materials and New Production Technologies** (Grant Agreement NMP2-SE-2009-228893) with a total cost of €5,305,678 and an EC contribution of €3,729,043. This is a 48 month project which started on 01-10-09.





A-FOOTPRINT Newsletter 5



which will then be scaled and refined, with a range of realworld foot and ankle pathological models added over time. The **ABT** team working with clinical and SME partners are also undertaking work to validate the models, and to assess the feasibility of driving the models using plantar pressure distribution data. Output data from the ABT modelling approach will also be used in by **MAS** in the MADYMO platform. This is a multi-body, finite-element software package will be used to predict motions and forces as well as tissue stresses and strains. This approach will be used to optimise or 'fine tune' the prototype orthotic designs developed by **Materialise** (**MAT**) and the **University of Newcastle upon Tyne** (**UNEW**).

The research work being undertaken in these two work packages is progressing very well. The UK clinical centre has gained ethical approval for subject and patient involvement and the first data sets have been collected and transferred to the relevant partners. The Dutch clinical centre will collect their first data sets once ethical approval has been granted. **ABT** have started the development of

Fig 3. 3D leg and foot reconstruction from CT scans

the generic foot model building in the individual bones, major ligaments and muscles from the anthropometric and CT data provided (Figure 4).

Over the following weeks and months the *A*-FOOTPRINT Consortium partners will start to gain new insights on the complex movement patterns and forces experienced in a highly complex structure that is the foot. Moreover, how disease processes and injuries determinentally alter function is urgently required to inform the design and manufacture of highly personalised orthotic devices.

A-FOOTPRINT CONTACT

Michelle Connolly *EU Project Executive Administrator* Tel: +44 (0)141 331 8956 Email: michelle.connolly@gcu.ac.uk




Fig 4. Preliminary foot model in the AnyBody Modelling platform



A-FOOTPRINT is a thirteen partner Consortium comprising six SMEs, five HEIs, one large-industry, and one research organisation partners from seven EU member states. The project is funded under **FP7 Cooperation Work Programme: Theme 4 - Nanosciences**, **Nanotechnologies**, **Materials and New Production Technologies** (Grant Agreement NMP2-SE-2009-228893) with a total cost of \in 5,305,678 and an EC contribution of \in 3,729,043. This is a 48 month project which started on 01-10-09.

