

**1 PATIENT** **2 PRESCRIBER** **3 ORDER INFO**

PATIENT NAME \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
DATE D / M / Y SEX M / F WEIGHT \_\_\_\_\_  
FOOTWEAR TYPE \_\_\_\_\_ SHOE SIZE \_\_\_\_\_

PRESCRIBER \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
PHONE \_\_\_\_\_  
EMAIL \_\_\_\_\_

DATE D / M / Y  
WO # \_\_\_\_\_  
 POST TO PATIENT  
REPEAT ORDER No  
 DUPLICATE (no change)  
 CHANGE (as marked)

**4 DEVICE TYPE**

- BRACE**
- RICHIE BRACE (standard brace)
    - Full Articulation
    - Permanently Fixed
  - RICHIE DYNAMIC ASSIST
    - Tamarac Hinge
  - RICHIE SOCCER (shin guard)
  - LITTLE RICHIE (1-12 years)
  - RICHIE GAUNTLET
    - 7" Most versatile height
    - 9" For maximum rigidity and control
  - RICHIE CALIFORNIA
    - 7" Most versatile height cast must measure 7"
    - 9" For maximum rigidity and control cast must measure 9"
  - RICHIE SOLID  
(Cast must capture full calf circumference)

- FOOT**
- LEFT  RIGHT  PAIR
  - Separate cast & separate order form enclosed for orthotic of non-braced foot.

**5 OPTIONAL SPECS**

- TOP COVER**
- MATERIAL**
- Puff  1.5mm  3mm  
Plastazote  1.5mm  3mm
- COMBINE MATERIAL WITH**
- Poron  1.5mm  3mm
- LENGTH**
- mets  sulcus  toes
- BOTTOM COVER SELECTION**
- MATERIAL**
- Agoflex
  - Vinyl
  - Puff 1.5mm
  - Nyplex 1.5mm
  - Microsuede
  - WSL Carbon

**6 POSTING, ADDITIONS & MODS**

- POSTING**
- POST TO REARFOOT**
- Calcaneal Vertical  Neutral (as they sit)
- POST TO FOREFOOT**
- extrinsic  intrinsic L \_\_\_° VR/VLG R \_\_\_° VR/VLG  
 to sulcus  to toes
- Skive L \_\_\_mm R \_\_\_mm
- CAST DRESSING**
- Skin  Type 1  Min  Moderate  Max

- SHELL MODIFICATIONS**
- L/R Heel Cup Depth  15mm  25mm  35mm
  - L/R 1st Met Cut-Out  w/ support post
  - L/R 1st Ray Cut-Out
  - L/R Rigid 1st Extension  sulcus  toes
  - L/R Medial Flange  standard  50%
  - L/R Lateral Flange  standard  50%
  - L/R Narrow Shell Grind
  - L/R Wide Shell Grind
  - L/R Fascial Accommodation
  - L/R Heel Spur Accommodation
  - L/R Heel Hole  w/gel plug  w/ poron plug
  - L/R Lateral Arch Suspender (not with lateral flange)
  - L/R Medial Arch Suspender (not with medial flange)
  - L/R Sweet Spot

- ADDITIONS**
- L/R Heel Lift \_\_\_ mm  EVA  poron  loose
  - L/R Heel Pad  1.5mm  3mm
  - L/R Morton's Extension
  - L/R Reverse Morton's Ext.  cushion cork  poron  EVA
  - L/R Functional Hallux Limitus Accommodation
  - Lesion Accommodation L [5] [4] [3] [2] [1] R [1] [2] [3] [4] [5]
  - L/R Horseshoe Spur
  - L/R Metatarsal Pad  Standard  50%
  - L/R Metatarsal Bar  Standard  50%
  - L/R Metatarsal Raise
  - L/R Poron Arch Pad  1.5mm  3mm
  - L/R Neuroma Pad
  - L/R Cuboid Pad

**EFFECTIVE PATHOLOGY-SPECIFIC PRESCRIPTIONS**

- Posterior Tibial Tendon Dysfunction**
- Original Richie Brace
  - Full articulation pivot
  - 4mm/6mm Medial skive
  - 35mm Heel cup
  - Medial flange optional & Sweet spot optional
- Lateral Ankle Instability**
- Original Richie Brace
  - Full articulation pivot
  - 2°/4° Valgus forefoot post to sulcus
  - 35mm Heel cup
- DJD of the Ankle/Subtalar Joint**
- Original Richie Brace
  - Permanently fixed pivot
  - Medial flange
  - 35mm Heel cup
- Dropfoot Deformity**
- Original Richie Brace / Dynamic
  - Permanently fixed or Tamarac pivot
  - 15mm/25mm Heel cup
- Note: The Richie Brace is not indicated for patients with total dorsiflexion paralysis, an unstable knee or if patient weighs over 250lbs.
- Charcot Deformity**
- Original Richie Brace
  - Permanently fixed pivot
  - Poron & Plastazote cover
  - 15mm/25mm Heel cup
- Note: The Richie Brace is not indicated for patients with severe "rocker bottom" deformity
- Caution:** The Richie Brace is not recommended if equinus is a primary deformity or for high risk diabetic patients

**7 SPECIAL INSTRUCTIONS**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



