

## 1 PATIENT

PATIENT NAME \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
D.O.B. D / M / Y SEX M / F WEIGHT \_\_\_\_\_  
FOOTWEAR TYPE \_\_\_\_\_ SHOE SIZE \_\_\_\_\_

## 2 PRESCRIBER

PRESCRIBER \_\_\_\_\_  
ADDRESS \_\_\_\_\_  
PHONE \_\_\_\_\_  
EMAIL \_\_\_\_\_

## 3 ORDER INFO

DATE D / M / Y  
WO # \_\_\_\_\_  
 RUSH  
 POST TO PATIENT  
REPEAT ORDER # \_\_\_\_\_  
 DUPLICATE (no change)  
 CHANGE (as marked)

## 4 DEVICE TYPE

**D** Diabetic **RA** Rheumatoid Arthritic  
Note: some devices are suitable for both Diabetic and Arthritic patients.

### ACCOMMODATIVE - Rubberflex Shell

CONTROL LEVEL 2-3

- D** STANDARD TRIDENSITY  
3.0mm P-Cell, 3.0mm Poron, Rubberflex & 50 durometer Fish Foam, 1.5mm Cushion Cork bottom cover.
- D**  **RA** SOFT MOULD  
1.5mm Puff, 1.5mm Poron, Rubberflex & Thermocork, Poron Fill, 1.5mm Nyplex bottom cover.
- RA** STANDARD TRIDENSITY  
3.0mm Bamboolon, Rubberflex & 50 durometer Fisher Foam, 1.5mm Nyplex bottom cover.
- D** CHARCOT TRIDENSITY  
3.0mm P-Cell, 3.0mm Poron, Rubberflex & Thermocork, 35 durometer Cloud EVA, 1.5mm Cushion Cork bottom cover.

### HYBRID DEVICES - Polypropylene shells [semi-flexible]

CONTROL LEVEL 3

- D** SEMI-FLEXIBLE MOULD  
3.0mm P-Cell, 3.0mm Poron, Polypropylene, Poron Fill, 1.5mm Cushion Cork bottom cover.
- RA** SEMI-FLEXIBLE MOULD  
3.0mm Bamboolon, Polypropylene, Poron Fill, 1.5mm Nyplex bottom cover.
- D**  **RA** SEMI-FLEX LOW PROFILE MOULD  
3.0mm Bamboolon, Polypropylene, Poron Fill, Microsuede bottom cover.

### FUNCTIONAL DEVICES - Polypropylene Shell [semi-rigid]

CONTROL LEVEL 4

- D** SEMI-RIGID MOULD  
3.0mm P-Cell, 3.0mm Poron, Polypropylene, Poron Fill, 1.5mm Cushion Cork bottom cover.
- RA** SEMI-RIGID MOULD  
3.0mm Bamboolon, Polypropylene, Poron Fill, 1.5mm Nyplex bottom cover.
- D**  **RA** SEMI-RIGID LOW PROFILE MOULD  
3.0mm Bamboolon, Polypropylene, Poron Fill, Microsuede bottom cover.

## 7 SPECIAL INSTRUCTIONS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_



## 5 OPTIONAL SPECS

### ULTRA PROTECTIVE TOP COVER

#### PREMIUM CUSHIONING FOR HIGH RISK FEET.

- Puff  1.5mm  3mm  
Plastazote  1.5mm  3mm

### FOREFOOT EXTENSION


#### MATERIAL

- Poron  1.5mm  3mm  
Nyplex  1.5mm  3mm  
Cushion Cork  1.5mm  3mm  
Puff  1.5mm  3mm

#### LENGTH

- sulcus  toes

### LESION ACCOMMODATION



### BOTTOM COVER SELECTION

#### MATERIAL

- Agoflex
- Vinyl
- Puff 1.5mm
- Nyplex 1.5mm
- Microsuede
- WSL Carbon

## 6 POSTING, ADDITIONS & MODS

### POSTING

#### POST TO

- Lab evaluation  Measurements from vertical  
 Calcaneal Vertical  Neutral (as they sit)

#### REARFOOT

- extrinsic  intrinsic L \_\_\_ ° VR/VLG R \_\_\_ ° VR/VLG

#### FOREFOOT

- extrinsic  intrinsic L \_\_\_ ° VR/VLG R \_\_\_ ° VR/VLG

- to sulcus  to toes

- Skive (needs deep heel cup) L \_\_\_ mm R \_\_\_ mm

- Inverted L \_\_\_ ° R \_\_\_ °

- Motion L \_\_\_ ° R \_\_\_ °

### CAST DRESSING

- Skin  Type 1  Min  Moderate  Max

### SHELL MODIFICATIONS

#### L/R Heel Cup Depth

- 14mm  16mm  18mm  20mm  22mm



- L/R 1st Met Cut-Out  w/ support post

- L/R 1st Ray Cut-Out

- L/R Rigid 1st Extension  sulcus  toes

- L/R Medial Flange  standard  50%

- L/R Lateral Flange  standard  50%

- L/R Narrow Shell Grind

- L/R Wide Shell Grind

- L/R Fascial Accommodation



- L/R Heel Spur Accommodation

- L/R Heel Hole  w/gel plug  w/ poron plug

### ADDITIONS

- L/R Heel Lift \_\_\_ mm  EVA  poron  loose



- L/R Heel Pad  1.5mm  3mm

- L/R Morton's Extension



- L/R Reverse Morton's Ext.  cushion cork  poron  EVA

- L/R Functional Hallux Limitus Accommodation



- L/R Horseshoe Spur

- L/R Metatarsal Pad  Standard  50%



- L/R Metatarsal Bar  Standard  50%

- L/R Metatarsal Raise



- L/R Poron Arch Pad  1.5mm  3mm

- L/R Neuroma Pad



- L/R Cuboid Pad

- L/R Amputation Accommodation  Partial  Full

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y

D/M/Y